

Product datasheet for **HP210601**

CHMP2B Human qPCR Primer Pair (NM_014043)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25978
Forward Sequence:	GGCTATAATCAGAGATCGAGCAG
Reverse Sequence:	CTCGTCTTCTGTTCCGTAGATG
ACCN:	<u>NM_014043</u> , <u>NM_014043.1</u> , <u>NM_014043.2</u> , <u>NM_014043.3</u> , <u>BC001553</u> , <u>BC001553.1</u> , <u>BI821848</u> , <u>BM828410</u> , <u>BM982148</u> , <u>NM_014043.4</u>
UniProt ID:	<u>Q9UQN3</u>
Synonyms:	ALS17; CHMP2.5; DMT1; FTDALS7; VPS2-2; VPS2B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)