

## Product datasheet for **HP210612**

### OIF (OGN) Human qPCR Primer Pair (NM\_014057)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 4969  |
| Forward Sequence: | CCATAATGCCCTGGAATCCGTG  |
| Reverse Sequence: | CAGGCGTATCTCTCAATGCGG   |
| ACCN:             | <a href="#">NM_014057</a> , <a href="#">NM_014057.1</a> , <a href="#">NM_014057.2</a> , <a href="#">NM_014057.3</a> , <a href="#">NM_014057.4</a> , <a href="#">BC095443</a> , <a href="#">BC095443.1</a> , <a href="#">BC037273</a> , <a href="#">BG772245</a> , <a href="#">BM679549</a> , <a href="#">BX646984</a> , <a href="#">NM_014057.5</a> |
| UniProt ID:       | <a href="#">P20774</a>  |
| Synonyms:         | OG; OIF; SLRR3A   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.        |
| Storage:          | Store at -20°C.   |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)