

Product datasheet for **HP210617**

NOB1 Human qPCR Primer Pair (NM_014062)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	28987
Forward Sequence:	GCACTCACATACCAGTTGGAAGC
Reverse Sequence:	GAGGTGTTTCTGGGTGCTGAATC
ACCN:	NM_014062 , NM_014062.1 , NM_014062.2 , BC064630 , BC064630.1 , BC000050 , BG386829 , BM793847 , BM827373 , NM_014062.3
UniProt ID:	Q9ULX3
Synonyms:	ART-4; MST158; MSTP158; NOB1P; PSMD8BP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)