

Product datasheet for **HP210618**

METTL11A (NTMT1) Human qPCR Primer Pair (NM_014064)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	28989
Forward Sequence:	CGGACTCTTACGACGTGATCTG
Reverse Sequence:	GCCATGTTGTCTTTGATGACGATG
ACCN:	<u>NM_014064</u> , <u>NM_014064.1</u> , <u>NM_014064.2</u> , <u>NM_014064.3</u> , <u>BC001396</u> , <u>BC001396.2</u> , <u>BC033234</u> , <u>BE392289</u> , <u>BF316784</u> , <u>BG468453</u> , <u>BI083998</u> , <u>BM009113</u> , <u>BM554138</u> , <u>BQ636193</u> , <u>BQ923165</u> , <u>NM_014064.4</u>
UniProt ID:	<u>Q9BV86</u>
Synonyms:	AD-003; C9orf32; HOMT1A; METTL11A; NRMT; NRMT1; NTM1A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)