

Product datasheet for **HP210659**

WDR91 Human qPCR Primer Pair (NM_014149)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29062
Forward Sequence:	AAACGGGACAGACTGCTCTTGC
Reverse Sequence:	CACGCAAGAGACAGGATTCTGG
ACCN:	NM_014149 , NM_014149.1 , NM_014149.2 , NM_014149.3 , BC017246 , BC017246.2 , BC009939 , BC063394 , BX640957 , NM_014149.4
UniProt ID:	A4D1P6
Synonyms:	HSPC049; SORF-1; SORF1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)