

Product datasheet for **HP210669**

MRPL18 Human qPCR Primer Pair (NM_014161)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29074
Forward Sequence:	TGGCACAGAGATGCTTAGAGGC
Reverse Sequence:	CACACCACCTTCTGTCATGGCA
ACCN:	NM_014161 , NM_014161.1 , NM_014161.2 , NM_014161.3 , NM_014161.4 , BC001623 , BC001623.1 , BF797811 , BG548185 , BM924328 , BU601065 , BU932963 , BU943173 , NM_014161.5
UniProt ID:	Q9H0U6
Synonyms:	HSPC071; L18mt; MRP-L18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)