

Product datasheet for **HP210677**

CHMP4A Human qPCR Primer Pair (NM_014169)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29082
Forward Sequence:	GGAGCAGAAGATTCAACAGGAGC
Reverse Sequence:	CGCTGAAACTCCAGGGTGGATA
ACCN:	<u>NM_014169</u> , <u>NM_014169.1</u> , <u>NM_014169.2</u> , <u>NM_014169.3</u> , <u>BC113533</u> , <u>BC010893</u> , <u>BC107699</u> , <u>BC113535</u> , <u>BC143834</u> , <u>BM041056</u>
UniProt ID:	<u>Q9BY43</u>
Synonyms:	C14orf123; CHMP4; CHMP4B; HSPC134; SHAX2; SNF7; SNF7-1; VPS32-1; VPS32A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)