

Product datasheet for **HP210718**

VAMP2 Human qPCR Primer Pair (NM_014232)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6844
Forward Sequence:	CTCCAAACCTCACCAGTAACAGG
Reverse Sequence:	AGCTCCGACAGCTTCTGGTCTC
ACCN:	NM_014232 , NM_014232.1 , NM_014232.2 , BC033870 , BC033870.1 , BC002737 , BC019608 , BU625107 , NM_014232.3
UniProt ID:	P63027
Synonyms:	NEDHAHM; SYB2; VAMP-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)