

Product datasheet for **HP210720**

UBL4A Human qPCR Primer Pair (NM_014235)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8266
Forward Sequence:	TGGTCTCCGAGAAGCTGAACGT
Reverse Sequence:	GACCACTAGGTTGAGCTTGGAG
ACCN:	<u>NM_014235</u> , <u>NM_014235.1</u> , <u>NM_014235.2</u> , <u>NM_014235.3</u> , <u>NM_014235.4</u> , <u>BC043346</u> , <u>BC020254</u> , <u>BC053589</u> , <u>BM711450</u> , <u>NM_014235.5</u>
UniProt ID:	<u>P11441</u>
Synonyms:	DX254E; DXS254E; G6PD; GDX; GET5; MDY2; TMA24; UBL4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)