

Product datasheet for **HP210724**

EIF2B2 Human qPCR Primer Pair (NM_014239)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8892
Forward Sequence:	CCAGTCCAACATCATTGAGGCG
Reverse Sequence:	ACTGTTCCGGGAGAAGCCAATGG
ACCN:	NM_014239 , NM_014239.1 , NM_014239.2 , NM_014239.3 , BC003165 , BC003165.1 , BC000494 , BC011750 , BI908906 , BX161473 , NM_014239.4
UniProt ID:	P49770
Synonyms:	EIF-2Bbeta; EIF2B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)