

Product datasheet for **HP210728**

ADAMTS2 Human qPCR Primer Pair (NM_014244)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9509
Forward Sequence:	TACAAGGACGCCTTCAGCCTCT
Reverse Sequence:	CCACTTTGCAGTGGCTGTTGTC
ACCN:	<u>NM_014244</u> , <u>NM_014244.1</u> , <u>NM_014244.2</u> , <u>NM_014244.3</u> , <u>NM_014244.4</u> , <u>BC098581</u> , <u>BC111459</u> , <u>BM512469</u> , <u>BU629258</u> , <u>BX479109</u> , <u>BX537606</u>
UniProt ID:	<u>O95450</u>
Synonyms:	ADAM-TS2; ADAMTS-2; ADAMTS-3; EDSDERMS; NPI; PC I-NP; PCI-NP; PCINP; PCPNI; PNPI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)