

Product datasheet for **HP210731**

RAPGEF2 Human qPCR Primer Pair (NM_014247)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9693
Forward Sequence:	CTCGGATCAGTATCTTGCCACAG
Reverse Sequence:	AGGTTCCACTGACAGGCAATGC
ACCN:	<u>NM_014247</u> , <u>NM_014247.1</u> , <u>NM_014247.2</u> , <u>BC110355</u> , <u>BC117321</u> , <u>BP386101</u> , <u>BU620827</u>
UniProt ID:	<u>Q9Y4G8</u>
Synonyms:	CNrasGEF; NRAPGEP; nRap GEP; PDZ-GEF1; PDZGEF1; RA-GEF; RA-GEF-1; RAGEF; Rap-GEP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)