

Product datasheet for **HP210739**

B3GNT3 Human qPCR Primer Pair (NM_014256)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10331
Forward Sequence:	AGGCACAGACTCACGGAGACAT
Reverse Sequence:	GTTGAGCACGAAGCTGGCGTTG
ACCN:	<u>NM_014256</u> , <u>NM_014256.1</u> , <u>NM_014256.2</u> , <u>NM_014256.3</u> , <u>BC067423</u> , <u>BC067423.1</u> , <u>BC066876</u> , <u>NM_014256.4</u>
UniProt ID:	<u>Q9Y2A9</u>
Synonyms:	B3GAL-T8; B3GN-T3; B3GNT-3; beta3Gn-T3; HP10328; TMEM3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)