

## Product datasheet for **HP210740**

### CD299 (CLEC4M) Human qPCR Primer Pair (NM\_014257)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10332
Forward Sequence:	GAGTAACCGCTTCTCCTGGATG
Reverse Sequence:	CGCACAGTCTTCATTCCCGCTA
ACCN:	<a href="#">BC038851</a> , <a href="#">NM_014257</a> , <a href="#">NM_014257.1</a> , <a href="#">NM_014257.2</a> , <a href="#">NM_014257.3</a> , <a href="#">NM_014257.4</a> , <a href="#">BC038851.1</a> , <a href="#">BC110614</a> , <a href="#">BQ008265</a> , <a href="#">BX102225</a> , <a href="#">NM_014257.5</a>
UniProt ID:	<a href="#">Q9H2X3</a>
Synonyms:	CD209L; CD299; DC-SIGN2; DC-SIGNR; DCSIGNR; HP10347; L-SIGN; LSIGN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)