

Product datasheet for **HP210741**

SCP2 (SYCP2) Human qPCR Primer Pair (NM_014258)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10388
Forward Sequence:	ACCAGATTCACAGGCAGCGGAA
Reverse Sequence:	GGTGCCAACATGCCCATTTGC
ACCN:	<u>NM_014258</u> , <u>NM_014258.1</u> , <u>NM_014258.2</u> , <u>NM_014258.3</u> , <u>BC040566</u> , <u>BC132864</u> , <u>BC132870</u> , <u>BC146798</u> , <u>BQ575787</u> , <u>BX648091</u> , <u>BX649108</u> , <u>NM_014258.4</u>
UniProt ID:	<u>Q9BX26</u>
Synonyms:	SCP-2; SCP2; SPGF1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)