

## Product datasheet for **HP210758**

### HSPA4L Human qPCR Primer Pair (NM\_014278)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22824
Forward Sequence:	CACAGTGATGCTCCAATGGAGAC
Reverse Sequence:	CCTCTTCTGGAGTGTGTTTCAGC
ACCN:	<a href="#">NM_014278</a> , <a href="#">NM_014278.1</a> , <a href="#">NM_014278.2</a> , <a href="#">NM_014278.3</a> , <a href="#">BC040560</a> , <a href="#">BC040560.1</a> , <a href="#">BM679019</a> , <a href="#">NM_014278.4</a>
UniProt ID:	<a href="#">O95757</a>
Synonyms:	APG-1; APG1; HSPH3; Osp94
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)