

Product datasheet for **HP210761**

HABP4 Human qPCR Primer Pair (NM_014282)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22927
Forward Sequence:	AAGAGCTGAGCGGAGATCCTAC
Reverse Sequence:	TCCTCTCAACGGTCTGTCTCGA
ACCN:	<u>NM_014282</u> , <u>NM_014282.1</u> , <u>NM_014282.2</u> , <u>NM_014282.3</u> , <u>BC018788</u> , <u>BI757571</u> , <u>BQ962055</u> , <u>NM_014282.4</u>
UniProt ID:	<u>Q5JVS0</u>
Synonyms:	IHABP-4; IHABP4; Ki-1/57; SERBP1L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)