

Product datasheet for **HP210763**

RRP4 (EXOSC2) Human qPCR Primer Pair (NM_014285)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23404
Forward Sequence:	CTCGCTGGTAACTCAGAGGATG
Reverse Sequence:	CTCCTGTTCCAAAAGCCTCTGG
ACCN:	NM_014285 , NM_014285.1 , NM_014285.2 , NM_014285.3 , NM_014285.4 , NM_014285.5 , NM_014285.6 , BC000747 , BC000747.2 , BC128445 , BU934039 , NM_014285.7
UniProt ID:	Q13868
Synonyms:	hRrp4p; p7; RRP4; Rrp4p; SHRF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)