

Product datasheet for **HP210766**

ITGB3BP Human qPCR Primer Pair (NM_014288)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23421
Forward Sequence:	CCACAAGTTCTGAAGAGCAAAAGC
Reverse Sequence:	CTCAACTTTTGATAGCAACATCATG
ACCN:	NM_014288 , NM_014288.1 , NM_014288.2 , NM_014288.3 , NM_014288.4 , BC009929 , BC009929.2 , BC014385 , BG548447 , BU633831 , BU853146 , NM_014288.5
UniProt ID:	Q13352
Synonyms:	CENP-R; CENPR; HSU37139; NRIF3; TAP20
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)