

## Product datasheet for **HP210770**

### CBX6 Human qPCR Primer Pair (NM\_014292)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23466
Forward Sequence:	AGAGGGAGCGTGAGCTGTATGG
Reverse Sequence:	GCTCGGCTTGACAGAGAAATGC
ACCN:	<a href="#">NM_014292</a> , <a href="#">NM_014292.1</a> , <a href="#">NM_014292.2</a> , <a href="#">NM_014292.3</a> , <a href="#">NM_014292.4</a> , <a href="#">BC012111</a> , <a href="#">BC012111.1</a> , <a href="#">BC064900</a> , <a href="#">NM_014292.5</a>
UniProt ID:	<a href="#">O95503</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)