

## Product datasheet for **HP210793**

### ORC6L (ORC6) Human qPCR Primer Pair (NM\_014321)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23594
Forward Sequence:	GAGAAGATTGGACAGCAGGTCG
Reverse Sequence:	GGTTTATGTGGCATCTCCTCTAC
ACCN:	<u><a href="#">NM_014321</a></u> , <u><a href="#">NM_014321.1</a></u> , <u><a href="#">NM_014321.2</a></u> , <u><a href="#">NM_014321.3</a></u> , <u><a href="#">BC063565</a></u> , <u><a href="#">BC063565.1</a></u> , <u><a href="#">BC039032</a></u> , <u><a href="#">BC050725</a></u> , <u><a href="#">BG822640</a></u> , <u><a href="#">NM_014321.4</a></u>
UniProt ID:	<u><a href="#">Q9Y5N6</a></u>
Synonyms:	ORC6L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)