

Product datasheet for **HP210803**

EID1 Human qPCR Primer Pair (NM_014335)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23741
Forward Sequence:	CCGCTCTTGAAGAAGCCGACAA
Reverse Sequence:	GGTCAGACGATTGACAACCATCA
ACCN:	<u>NM_014335</u> , <u>NM_014335.1</u> , <u>NM_014335.2</u> , <u>BC114944</u> , <u>BC007944</u> , <u>BC024151</u> , <u>BC114946</u> , <u>BI550673</u> , <u>NM_014335.3</u>
UniProt ID:	<u>Q9Y6B2</u>
Synonyms:	C15orf3; CRI1; EID-1; IRO45620; PNAS-22; PTD014; RBP21
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)