

Product datasheet for **HP210839**

DNAJC2 Human qPCR Primer Pair (NM_014377)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27000
Forward Sequence:	GCTAAACGGAAGGAGCAAGAAGC
Reverse Sequence:	CCTTCTTTGCCAGCAATGCTTGC
ACCN:	NM_014377 , NM_014377.1 , BC000859 , BC012376 , BC032854 , BC046351 , BC056682 , BC062703 , BC125057 , BC139751 , BC160045 , NM_014377.3
UniProt ID:	Q99543
Synonyms:	MPHOSPH11; MPP11; ZRF1; ZUO1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)