

Product datasheet for **HP210846**

UTP25 Human qPCR Primer Pair (NM_014388)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27042
Forward Sequence:	ACCTGATGCAGAACTGGGAGCA
Reverse Sequence:	ATTGAGGCTCCACATCCGCACT
ACCN:	NM_014388 , NM_014388.1 , NM_014388.2 , NM_014388.3 , NM_014388.4 , NM_014388.5 , NM_014388.6 , BC022964 , BC022964.1 , BX345520 , BX648514
UniProt ID:	Q68CQ4
Synonyms:	C1orf107; DEF; DIEXF; DJ434O14.5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)