

Product datasheet for **HP210858**

UQCRQ Human qPCR Primer Pair (NM_014402)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27089
Forward Sequence:	GCCTATCCGCACGTCTTCACTA
Reverse Sequence:	GGATCTCTCGAACTCTTCAGTCC
ACCN:	<u>NM_014402</u> , <u>NM_014402.1</u> , <u>NM_014402.2</u> , <u>NM_014402.3</u> , <u>NM_014402.4</u> , <u>BC090048</u> , <u>BC001390</u> , <u>BC064949</u> , <u>BG678118</u> , <u>BU190580</u> , <u>BU953465</u> , <u>NM_014402.5</u>
UniProt ID:	<u>O14949</u>
Synonyms:	MC3DN4; QCR8; QP-C; QPC; UQCR7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)