

Product datasheet for **HP210863**

KCNMB3 Human qPCR Primer Pair (NM_014407)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27094
Forward Sequence:	TCAGCCATCCAGGTCAGAAAGC
Reverse Sequence:	TCTATCTTGGTGGCACTTAGGTG
ACCN:	<u>NM_014407</u> , <u>NM_014407.1</u> , <u>NM_014407.2</u> , <u>NM_014407.3</u> , <u>BC082272</u> , <u>BC113851</u> , <u>BC152860</u> , <u>BM745010</u>
UniProt ID:	<u>Q9NPA1</u>
Synonyms:	BKBETA3; HBETA3; K(VCA)BETA-3; KCNMB2; KCNMBL; SLO-BETA-3; SLOBETA3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)