

## Product datasheet for **HP210866**

### SND1-IT1 Human qPCR Primer Pair (NM\_014411)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27099
Forward Sequence:	CAGCACATCAGTCATCGCCCTT
Reverse Sequence:	TGCTGTGGGAAGTATGGAAGG
ACCN:	<a href="#">NM_014411</a> , <a href="#">NM_014411.1</a> , <a href="#">BC093748</a> , <a href="#">BC112141</a>
Synonyms:	brain and nasopharyngeal carcinoma susceptibility protein; chromosome 7 open reading frame 54; MGC138346; NAG8; NAG8, NSG-X, MGC138346; nasopharyngeal carcinoma associated gene protein-8; NSG-X
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)