

Product datasheet for **HP210882**

IL20 Receptor alpha (IL20RA) Human qPCR Primer Pair (NM_014432)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	53832
Forward Sequence:	GCTGAAAGTGGACGGTTCTATCC
Reverse Sequence:	TTCCAATTCTCTGGAGCTGTCAG
ACCN:	<u>NM_014432</u> , <u>NM_014432.1</u> , <u>NM_014432.2</u> , <u>NM_014432.3</u> , <u>BC113574</u> , <u>BC113602</u> , <u>BM713319</u>
UniProt ID:	<u>Q9UHF4</u>
Synonyms:	CRF2-8; IL-20R-alpha; IL-20R1; IL-20RA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)