

## Product datasheet for **HP210887**

### IL37 Human qPCR Primer Pair (NM\_014439)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27178
Forward Sequence:	TCAGCCTCTGCGGAGAAAGGAA
Reverse Sequence:	GCGTGCTGATTCCTTTTGGGCA
ACCN:	<u><a href="#">NM_014439</a></u> , <u><a href="#">NM_014439.1</a></u> , <u><a href="#">NM_014439.2</a></u> , <u><a href="#">NM_014439.3</a></u> , <u><a href="#">BC020637</a></u>
UniProt ID:	<u><a href="#">Q9NZH6</a></u>
Synonyms:	FIL1; FIL1(ZETA); FIL1Z; IL-1F7; IL-1H; IL-1H4; IL-1RP1; IL-23; IL-37; IL1F7; IL1H4; IL1RP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)