

## Product datasheet for **HP210892**

### GCP4 (TUBGCP4) Human qPCR Primer Pair (NM\_014444)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 27229   |
| Forward Sequence: | CTTGGTGGACTTTGAACAGGTGG   |
| Reverse Sequence: | GCCTGAAACAGTTCTCCACGTC  |
| ACCN:             | <u><a href="#">BC000966</a></u> , <u><a href="#">BC009870</a></u> , <u><a href="#">NM_014444</a></u> , <u><a href="#">NM_014444.1</a></u> , <u><a href="#">NM_014444.2</a></u> , <u><a href="#">NM_014444.3</a></u> , <u><a href="#">NM_014444.4</a></u> , <u><a href="#">BC000966.2</a></u> , <u><a href="#">BC009870.2</a></u> , <u><a href="#">BC012801</a></u> , <u><a href="#">BC012801.2</a></u> , <u><a href="#">NM_014444.5</a></u> |
| UniProt ID:       | <u><a href="#">Q9UGJ1</a></u>   |
| Synonyms:         | 76P; GCP-4; GCP4; Grip76; MCCRP3  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.  |
| Storage:          | Store at -20°C.   |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)