

## Product datasheet for **HP210911**

### SRPX2 Human qPCR Primer Pair (NM\_014467)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27286
Forward Sequence:	CGTTACTACTGCCTATGACCGAG
Reverse Sequence:	GGTATTCACAGGTGGCACCATAG
ACCN:	<a href="#">NM_014467</a> , <a href="#">NM_014467.1</a> , <a href="#">NM_014467.2</a> , <a href="#">BC020733</a> , <a href="#">BC020733.1</a> , <a href="#">BP279501</a> , <a href="#">BU625416</a> , <a href="#">BX339682</a> , <a href="#">NM_014467.3</a>
UniProt ID:	<a href="#">O60687</a>
Synonyms:	BPP; CBPS; PMGX; RESDX; SRPUL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)