

Product datasheet for **HP210917**

DIMT1L (DIMT1) Human qPCR Primer Pair (NM_014473)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27292
Forward Sequence:	TGTTCAACACGGGGATTGGGCA
Reverse Sequence:	CAGGTCCAACCTCCAGCACTAC
ACCN:	<u>NM_014473</u> , <u>NM_014473.1</u> , <u>NM_014473.2</u> , <u>NM_014473.3</u> , <u>BC010874</u> , <u>BC002841</u> , <u>BC017885</u> , <u>BM749030</u>
UniProt ID:	<u>Q9UNQ2</u>
Synonyms:	DIM1; DIMT1L; HSA9761; HUSSY5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)