

Product datasheet for **HP210944**

APOBEC3C Human qPCR Primer Pair (NM_014508)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27350
Forward Sequence:	CGGAACGAAACTTGGCTGTGCT
Reverse Sequence:	CCAAGAGAGGAAGCACCTTTCTG
ACCN:	<u>NM_014508</u> , <u>NM_014508.1</u> , <u>NM_014508.2</u> , <u>BC011739</u> , <u>BC011739.2</u> , <u>BC021080</u>
UniProt ID:	<u>Q9NRW3</u>
Synonyms:	A3C; APOBEC1L; ARDC2; ARDC4; ARP5; bK150C2.3; PBI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)