

Product datasheet for **HP210962**

Grainyhead like protein 1 homolog (GRHL1) Human qPCR Primer Pair (NM_014552)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29841
Forward Sequence:	AACCGCAGCAACAAGCCTGTGC
Reverse Sequence:	TGAGAGGGAAGCAGTGCCACTT
ACCN:	NM_014552 , NM_014552.2 , BC067519 , BC067519.1 , BC063299 , BC067520 , BC067521 , BC098456 , BE537504 , BF573452
UniProt ID:	Q9NZI5
Synonyms:	grainyhead-like 1 (Drosophila); LBP-32; LBP32; LBP protein 32; leader-binding protein 32; mammalian grainyhead; MGR; TF2P2L2; transcription factor CP2-like 2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)