

Product datasheet for **HP210970**

OR1D5 Human qPCR Primer Pair (NM_014566)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8386
Forward Sequence:	TGGCTGCTTCATCTTCCTCACC
Reverse Sequence:	CCAAATGGGAGGCACAGGTAGA
ACCN:	NM_014566 , NM_014566.1 , BC106876 , BC152939
UniProt ID:	P58170
Synonyms:	C17orf2; OR1D4; OR17-2; OR17-30; OR17-31
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)