

Product datasheet for **HP210982**

SLC39A2 Human qPCR Primer Pair (NM_014579)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29986
Forward Sequence:	GGTTCAGATTGATGCAGCCAG
Reverse Sequence:	TTTCCTCCAGGGCTTCAGCAGT
ACCN:	<u>NM_014579</u> , <u>NM_014579.1</u> , <u>NM_014579.2</u> , <u>NM_014579.3</u> , <u>BC096723</u> , <u>BC098254</u> , <u>BC098274</u> , <u>BC110057</u> , <u>BP359281</u>
UniProt ID:	<u>Q9NP94</u>
Synonyms:	6A1; ETI-1; ZIP-2; ZIP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)