

## Product datasheet for **HP210988**

### SLC40A1 Human qPCR Primer Pair (NM\_014585)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	30061
Forward Sequence:	GAGACAAGTCCTGAATCTGTGCC
Reverse Sequence:	TTCTTGCAGCAACTGTGCACAG
ACCN:	<u><a href="#">NM_014585</a></u> , <u><a href="#">NM_014585.1</a></u> , <u><a href="#">NM_014585.2</a></u> , <u><a href="#">NM_014585.3</a></u> , <u><a href="#">NM_014585.4</a></u> , <u><a href="#">BC037733</a></u> , <u><a href="#">BC035893</a></u> , <u><a href="#">BX537503</a></u> , <u><a href="#">NM_014585.6</a></u>
UniProt ID:	<u><a href="#">Q9NP59</a></u>
Synonyms:	FPN1; HFE4; IREG1; MST079; MSTP079; MTP1; SLC11A3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)