

## Product datasheet for **HP210993**

### ERVWE1 (ERVW-1) Human qPCR Primer Pair (NM\_014590)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	30816
Forward Sequence:	CGCAACTGCTATCACTCTGCCA
Reverse Sequence:	CAGACAGTGACTCCAAGTCTC
ACCN:	<u><a href="#">NM_014590</a></u> , <u><a href="#">NM_014590.1</a></u> , <u><a href="#">NM_014590.2</a></u> , <u><a href="#">NM_014590.3</a></u> , <u><a href="#">BC137381</a></u> , <u><a href="#">BG012022</a></u> , <u><a href="#">BX365066</a></u> , <u><a href="#">BX391741</a></u>
UniProt ID:	<u><a href="#">Q9UQF0</a></u>
Synonyms:	ENV; ENVW; ERVWE1; HERV-7q; HERV-W-ENV; HERV7Q; HERVW; HERVWENV
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)