

Product datasheet for **HP211032**

PLCH2 Human qPCR Primer Pair (NM_014638)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9651
Forward Sequence:	TGGACCTGTCATCAGTGAGCAG
Reverse Sequence:	TCATCAGACACCTCGCCTTCCT
ACCN:	NM_014638 , NM_014638.1 , NM_014638.2 , BC019679 , BC037356 , BC043358 , BC044630 , BC050037 , BC128207 , BC128208 , BC144620 , BM971622 , BX117873
UniProt ID:	O75038
Synonyms:	PLC-eta2; PLC-L4; PLCeta2; PLCL4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)