

Product datasheet for **HP211033**

TTC37 Human qPCR Primer Pair (NM_014639)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9652
Forward Sequence:	GACACTGATGCTGAATCTGGAGC
Reverse Sequence:	CCATTTTGCCGTTCCAGCACTTG
ACCN:	<u>NM_014639</u> , <u>NM_014639.1</u> , <u>NM_014639.2</u> , <u>NM_014639.3</u> , <u>BC056893</u> , <u>BC056893.1</u> , <u>BC009767</u> , <u>BC015163</u> , <u>BC030966</u> , <u>NM_014639.4</u>
UniProt ID:	<u>Q6PGP7</u>
Synonyms:	KIAA0372; Ski3; THES
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)