

Product datasheet for **HP211043**

Importin 13 (IPO13) Human qPCR Primer Pair (NM_014652)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9670
Forward Sequence:	CTGTCAGCTCTTCAAGTGGAGG
Reverse Sequence:	TCCAGTGTGGTGAAGAGTTGG
ACCN:	NM_014652 , NM_014652.1 , NM_014652.2 , NM_014652.3 , BC008194 , BC008194.1 , BC002608 , BX116797
UniProt ID:	O94829
Synonyms:	IMP13; KAP13; LGL2; RANBP13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)