

Product datasheet for **HP211053**

CLINT1 Human qPCR Primer Pair (NM_014666)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9685
Forward Sequence:	CGGAGGATTTGCTGACTTTGGC
Reverse Sequence:	AACTCGCCACTGGAAGCAACAG
ACCN:	NM_014666 , NM_014666.1 , NM_014666.2 , NM_014666.3 , BC004467 , BC013091 , BC030264 , BI461076 , BK000414 , BP246338 , NM_014666.4
UniProt ID:	Q14677
Synonyms:	CLINT; ENTH; EPN4; EPNR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)