

Product datasheet for **HP211059**

MRPP3 (KIAA0391) Human qPCR Primer Pair (NM_014672)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9692
Forward Sequence:	GGAAGCACATGCTAAGACGGAG
Reverse Sequence:	GTGCAGTGTGGCATACAGAAGG
ACCN:	<u>NM_014672</u> , <u>NM_014672.1</u> , <u>NM_014672.2</u> , <u>NM_014672.3</u> , <u>BC032221</u> , <u>BC032221.1</u> , <u>BC001542</u> , <u>BC044580</u> , <u>BP373829</u> , <u>BQ063741</u> , <u>BX161394</u> , <u>BX161487</u>
UniProt ID:	<u>O15091</u>
Synonyms:	KIAA0391; MRPP3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)