

Product datasheet for **HP211064**

CEP57 Human qPCR Primer Pair (NM_014679)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9702
Forward Sequence:	GGAGTTGTCAGAAGTCTTACAGAC
Reverse Sequence:	CATCCTTCCCACTAATGCCTCC
ACCN:	BC009053 , NM_014679 , NM_014679.1 , NM_014679.2 , NM_014679.3 , NM_014679.4 , BC001233 , BC029385 , BC031048 , BC039711 , BG699097 , NM_014679.5
UniProt ID:	Q86XR8
Synonyms:	MVA2; PIG8; TSP57
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)