

Product datasheet for **HP211071**

Rubicon (RUBCN) Human qPCR Primer Pair (NM_014687)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9711
Forward Sequence:	CGACCTGGAAAAGGAGAATGCC
Reverse Sequence:	TCCTCTAGGCACTGGCTCATCA
ACCN:	BC033615 , NM_014687 , NM_014687.1 , NM_014687.2 , NM_014687.3 , BC033615.1 , BC014173 , BC160011 , BU501990 , BX647598 , NM_014687.4
UniProt ID:	Q92622
Synonyms:	KIAA0226; RUBICON; SCAR15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)