

Product datasheet for **HP211078**

CAPON (NOS1AP) Human qPCR Primer Pair (NM_014697)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9722
Forward Sequence:	AGAGACACCGCTGTCCACTCAC
Reverse Sequence:	AGCCAACGGTCCTTCAGCAAG
ACCN:	<u>NM_014697</u> , <u>NM_014697.1</u> , <u>NM_014697.2</u> , <u>BC143771</u> , <u>BC041382</u> , <u>BC111928</u> , <u>BC112295</u> , <u>BC118499</u> , <u>BC143780</u> , <u>NM_014697.3</u>
UniProt ID:	<u>O75052</u>
Synonyms:	6330408P19Rik; CAPON; NPHS22
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)