

Product datasheet for **HP211095**

ARHGAP32 Human qPCR Primer Pair (NM_014715)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9743
Forward Sequence:	GCAGTCAGAGTAAGGCTGTAGC
Reverse Sequence:	GGCAACATTTTTTCGGAGGCGGT
ACCN:	<u>NM_014715</u> , <u>NM_014715.1</u> , <u>NM_014715.2</u> , <u>NM_014715.3</u> , <u>BC104898</u> , <u>BC000277</u> , <u>BC023644</u> , <u>BC051236</u> , <u>BC054488</u> , <u>BC065371</u> , <u>BC113429</u> , <u>BQ319051</u> , <u>NM_014715.4</u>
UniProt ID:	<u>A7KAX9</u>
Synonyms:	GC-GAP; GRIT; p200RhoGAP; p250GAP; PX-RICS; RICS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)