

Product datasheet for **HP211123**

HURP (DLGAP5) Human qPCR Primer Pair (NM_014750)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9787
Forward Sequence:	CTCGATCAGCTACTCAAGCAGC
Reverse Sequence:	CAGGTCTTCCTTACTTGGCACC
ACCN:	NM_014750 , NM_014750.1 , NM_014750.2 , NM_014750.3 , NM_014750.4 , BC010658 , BC010658.1 , BC016276 , BM788630 , BT007344 , BX248255 , NM_014750.5
UniProt ID:	Q15398
Synonyms:	DLG7; HURP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)