

## Product datasheet for **HP211134**

### MRPL19 Human qPCR Primer Pair (NM\_014763)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9801
Forward Sequence:	CAGGAAGAGGACTTGGAGCTAC
Reverse Sequence:	GCTATCATCCAGCCGTTTCTCTA
ACCN:	<a href="#">NM_014763</a> , <a href="#">NM_014763.1</a> , <a href="#">NM_014763.2</a> , <a href="#">NM_014763.3</a> , <a href="#">BC030144</a> , <a href="#">BC030144.2</a> , <a href="#">BC021051</a> , <a href="#">BI551514</a> , <a href="#">BQ432684</a> , <a href="#">NM_014763.4</a>
UniProt ID:	<a href="#">P49406</a>
Synonyms:	L19mt; MRP-L15; MRP-L19; MRPL15; RLX1; RPML15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)