

Product datasheet for **HP211142**

DELE1 Human qPCR Primer Pair (NM_014773)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9812
Forward Sequence:	CTACCTGGATGAGCAGAGAGCT
Reverse Sequence:	CCAGATTCCTCTGCACACCAAG
ACCN:	<u>NM_014773</u> , <u>NM_014773.1</u> , <u>NM_014773.2</u> , <u>NM_014773.3</u> , <u>NM_014773.4</u> , <u>BC007855</u> , <u>BC007855.2</u> , <u>BC010452</u> , <u>BC011269</u> , <u>BC011926</u> , <u>BC035486</u>
UniProt ID:	<u>Q14154</u>
Synonyms:	DELE; DELE1(L); KIAA0141
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

